



50

50 Years of excellence
to microbiology

Ready to use Contact and Settle Plates



A wide range of QC Plates is available to meet the requirements of many industries and applications including

FOOD
DAIRY
BREWING
CLINICAL
COSMETICS
PHARMACEUTICAL
HEAVY INDUSTRY

QC Plates

Medical Wire & Equipment Co (Bath) Ltd

Corsham, Wiltshire, England SN13 9RT Tel: +44 (0) 1225 810361 Fax: +44 (0) 1225 810153 E-mail: info@mwe.co.uk website: www.mwe.co.uk

QC Plates

A range of ready to use contact and settle plates with a choice of medium to cover most applications.

- ◆ **Six-month shelf life !**
- ◆ Wide range of media to meet the requirements of many industries.
- ◆ Easy to use.
- ◆ Grid on base allows rapid and accurate analysis
- ◆ Individually sealed means no wastage – economical!



DipKit[®] TSA

Dip Kit[®] for Water and Oil Testing

Dip slide with Tryptone Soy Agar for total viable count determination for water, beverages, cutting oil, slurry, etc.

Dip Kits can also be supplied with CLED Agar or Maconkey Agar for more detailed analysis of microbial population (Please contact us for details)



Product Code	Product Description	Unit size
QC Plates		
MW600	Enterobacteriaceae (C.L.E.D.)	120
MW605	Yeasts & Moulds (mycomedia + chloramphenicol)	120
MW610	Hygiene Control (TSA + surfactant)	120
MW615	Coliform Selective Media	120
MW625	Listeria Selective Media (MOX)	120
MW630	Pseudomonas Selective Media	120
MW635	Total Viable Count (nutrient agar)	120
MW640	Bacillus Cereus Selective Media	120
Dip Kits[®]		
MW203	Dip Kit (TSA)	210

QC Plates